ISO 13485:2016

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OGM SPACER WITH ANTIBIOTIC KNEE, HIP, SHOULDER USER'S GUIDE

PRODUCT DESCRIPTION: Spacer products (OGM Hip Spacer with Antibiotic, OGM Knee Spacer with Antibiotic, OGM Shoulder Spacer with Antibiotic) are produced from OGM 3A Bone Cement with Antibiotic product. Mixed bone cement is filled in to the spacer mold by injection. In the spacer, are putted stainless steel to be improved resistance because of increased in hip and knee spacer length. Bone cement is made injection on this part and this part was covered completely. It occurs by direct casting of bone cement into the mold in the knee prosthesis.

COMPOSITION:

The powder component comprising polymethyl methacrylate, barium sulphate, benzoyl peroxide and gentamicin; Methyl Methacrylate is formed by mixing with the liquid component consisting of N,N dimethyl p-toluidine and Hydroquinone and pouring into the spacer mold by injection.

INTENDED USE AND INDICATIONS OF THE DEVICE:

It is a medical device that can be implanted temporarily. In chronic periprosthetic joint infections, systemic antibiotics and localized antibiotic-releasing temporary spacers are used following removal of the infected implant in the first stage. Spacers are useful in maintaining soft tissue blood pressure.

Gentamicin Properties in OGM Spacer: (i) Wide antibacterial spectrum, (ii) good bactericidal effect at low concentrations, (iii) low resistance development, (iv) low binding rate to proteins, (v) low allergic potential, (vi) ineffective or insignificant effect on bone cement-water mechanics, (vii) resistant to chemical and thermal factors, (viii) good solubility in water, (ix) good release from the spacer.

Advantages:

- Preformed shapes and ready to use in different sizes. Makes surgery quick.
- To create a joint cavity after removal of an infected prosthesis and to facilitate re-implantation by preventing contraction of muscles.
- Restoring joint functionality within the time required for complete drying of the infection, depending on specific clinical conditions, and providing the patient with a higher quality of life.
- Shortening the functional recovery time to provide significant fitness to the muscles holding the joint after the fixed prosthesis
- Devices have standardized and approved mechanical and pharmacological performances.

- Long-term antibiotic release.

CONTRAINDICATIONS OF THE DEVICE:

- The patient's condition is contraindicated for this two-stage arthroplasty procedure due to decreased immune response or other relevant systemic clinical conditions
- Lack of proper bone structure to achieve adequate property of prosthesis
- The procedure is not appropriate due to disorders of the muscle, nervous or vascular system of the patient
- Low bone quality causes replacement of the prosthesis and bone fractures
- Infected THR or TKR device cannot be removed
- Resistance to gentamicin of infection pathogens
- The patient is susceptible to gentamicin, aminoglycosides or PMMA bone filling
- Suspicion of a separate systemic or secondary infection or verification
- The patient does not have sufficient bone stock for insertion and fixation of the prosthesis
- Weight, age or activity level of the patient cause the surgeon to believe that he / she will experience early system disorder
- The patient does not want or cannot use protected weight bearing mobility (eg crutches, walking stick, walker, etc.) during implantation

SAFETY AND WARNINGS

Carefully read the instructions for use. Aseptic surgical techniques should be followed. For the safe application of OGM Spacer, surgeons must be familiar with product features, application limitations, and correct application techniques. Changing the structural identity of the device may damage the device (may cause cracks or fractures and may result in the loss of sterility, creating a risk of infection for the patient). Surgeons should ensure that the correct size device is properly installed at the implantation site. Remove the OGM Spacer not later than six months after implantation. The device has been tested for safe use for more than 6 months. Exceeding this period for a long time can lead to abrasion, residue formation and breakage.

Using OGM Spacer under conditions other than those specified does not benefit the patient and increases the risk of drug-resistant bacteria.

Do not resterilize and / or reuse the device. Designed on a single patient and for single use. The device should not be resterilized unless there is a risk of infection on the patient. Resterilization may alter the morphology of the device, the effectiveness or mechanical properties of antibiotics, and may pose serious risks to the patient's health. In any case, re-use of the device after removal should be prevented as it may cause contamination and worsen patient infection. During the removal, the spacer may be damaged or fill debris may remain on the device.

APPLICATION

After removal of the prosthesis and before insertion, the area should be washed with Ringer or physiological solution and the residues should be thoroughly cleaned.

In particular, any cement residues from the previous implantation should be removed.

To ensure optimal antibiotic release, it is necessary to prevent the insert from being washed with water solutions before and after implantation.

STORAGE

Store in a dark ,cool and dry place below 25 ° C.

Flammable - Keep away from sources of ignition. Humidity should be between 0-50.

SHELF LIFE / STERILITY

Shelf life is 3 years for sterile products.

Do not use the product if the expiry date indicated on the product package has passed.

Opened or damaged packages cannot be resterilized and must therefore be discarded. They were sterilized by gasification with spacer ethylene oxide.

PATIEN POPULATION

Suitable for all patient populations except for those who have a known or assumed hypersensitivity to the substances contained in bone cement and who are in pregnancy and lactation.

It is not used in patients who have an active or untreated infection in the bone area, due to non-gentamicin-inducing offenses.

OGM Spacer should not be used during the first trimester of pregnancy, and should not be used only in life-threatening situations for the remaining period. Its use in children should only be used when it is not possible to interfere with other procedures.






RESPONSIBILITY

The manufacturer shall not be liable for any indirect or penal sanctions or transactions arising from or associated with the re-use of bone cement.

INTENDED USER

Designed for use by healthcare professionals

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