2195

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OGM VTB STERILE RADIOPAQUE VERTEBROPLASTY **CEMENT USER GUIDE**

PRODUCT DESCRIPTION:

OGM VTB Sterile Radiopaque Vertebroplasty Bone Cement is injected into a partially collapsed vertebral body (fractured vertebra). The cement strengthens the bone, reduces fracture, alleviates pain, and is used to provide mechanical stability in cases of osteoporosis or tumor infiltration.

It is an acrylic compound containing powder made of methyl methacrylate polymer combined with benzoyl peroxide and liquid monomer methyl methacrylate activated with N,N-dimethyl-p-toluidine, capable of self-polymerization. When the monomer and polymer are mixed appropriately, N, N-dimethyl-p-toluidine activates the benzoyl peroxide and initiates polymerization of the monomer.

The polymerization process occurs progressively. The mixture is initially in liquid form, then becomes creamy, transforms into a dough-like consistency, and finally hardens. This reaction is exothermic.

COMPOSITION:

The package contains one bag of cement powder (polymer powder) and one colored glass ampoule (monomer liquid).

- **Powder component:** Contains polymethyl methacrylate (Acrylic Polymer), barium sulfate, and benzoyl peroxide.
- **Liquid component:** Contains methyl methacrylate, N,N-dimethyl-p-toluidine, and hydroquinone.

Vertebroplasty Bone Cement

Formula of Components	OGM VTB 12,5 Vertebroplasty Cement	OGM VTB 15 Vertebroplasty Cement	OGM VTB 20 Vertebroplasty Cement	OGM VTB 25 Vertebroplasty Cement
	REF 0308/VTB	REF 2023/VTB	REF 2071/VTB	REF 1406/VTB

<u>Liquid Component</u>	5 ml ± 0.25 ml	6 ml ± 0.3 ml	8 ml ± 0.4 ml	10 ml ± 0.5 ml
Methyl Methacrylate	98% w/w	98% w/w	98% w/w	98% w/w
N,Ndimethyl-p-toluidene	2% w/w	2% w/w	2% w/w	2% w/w
Hydroquinone	20 ppm	20 ppm	20 ppm	20 ppm
<u>Powder Component</u>	12.5 g ± 0.625 g	15 g ± 0.75 g	20 g ± 1 g	25 g ± 1.25 g
Polymethyl methacrylate	59.36% w/w	59.36% w/w	86.8% w/w	86.8% w/w
Benzoyl Peroxide	0.64% w/w	0.64% w/w	0.64% w/w	0.64% w/w
Barium Sulfate	40% w/w	40% w/w	40% w/w	40% w/w

INTENDED USE AND INDICATIONS:

OGM VTB Vertebroplasty is indicated for the following conditions:

- Compression fractures of the vertebral body
- Vertebral body tumors (metastases or myelomas)
- Symptomatic vertebral hemangiomas
- Osteoporosis
- Tumor infiltration
- It is designed to fill vertebral bodies with secondary structural changes in selected trauma cases, and its use is recommended in these situations.

CONTRAINDICATIONS:

- Use is contraindicated in patients with known allergy to any of the product components.
- Hemorrhagic diathesis and active infection are absolute contraindications for vertebroplasty bone cement.
- Vertebral body lesions with epidural extension are considered a relative contraindication due to the risk of spinal cord compression.

COMPLICATIONS:

- Superficial and deep infection
- Cement leakage into blood vessels
- Mobilization of bone fragments during injection
- Cement leakage into the perivertebral space
- Potential involvement of the spinal cord, nerve roots, peritoneal and intrathoracic structures, or retroperitoneal region

SIDE EFFECTS:

The following side effects and complications may occur after vertebroplasty procedures:

- Leakage of cement into the vertebral canal or adjacent vessels
- Superficial or deep infection
- Bleeding
- Increased back pain
- Fracture of the vertebra or adjacent ribs
- Numbness or tingling in the spine
- Paralysis
- Vasovagal reaction
- Cardiac rhythm changes
- Acute myocardial infarction
- Cardiorespiratory arrest
- Cerebrovascular events
- Pulmonary thromboembolism
- Mobilization of bone fragments during injection

- Cement leakage into the perivertebral space
- Potential involvement of the spinal cord, nerve roots, peritoneal and intrathoracic structures, or the retroperitoneal region

During application, modern cementing techniques should be followed to limit undesired side effects and ensure secure, long-term anchoring of the bone cement within the vertebral body.

SAFETY AND WARNINGS:

The monomer is a volatile liquid and potentially irritant; it must be handled with care:

- Avoid unnecessary contact with skin and mucous membranes. The use of double gloves is recommended.
- Avoid exposure to vapors generated during cement mixing.
- When opening the ampoule, break the neck carefully to prevent cuts or splashing. After opening, pour the contents quickly into the mixing container and add the polymer.
- Cement and mixing components should not be kept above 25 °C. Higher temperatures shorten the working time, while lower temperatures extend it.
- Air entry into the cement must be prevented. Otherwise, air bubbles may form within the dough-like mixture and at the bone interface. This must be avoided during mixing, transfer, and application.
- This procedure should only be performed by medical professionals experienced in the technique, in a healthcare facility equipped with all necessary equipment. This ensures correct application and allows prompt intervention in case of any undesired patient reactions.
- Increased application pressure may potentially damage the vertebral body and surrounding soft tissues.
- Viscosity increases progressively after mixing. After a certain period and as the cement hardens, an increase in viscosity is normal. This increase raises resistance to application and reduces cement pressure within the vertebra.

Pharmaceutical Warnings:

- Do not use products past the expiration date indicated on the label.
- Single-use product; intended for one patient only. Reuse may pose mechanical, physicochemical, and biological contamination risks.

- Do not use opened or damaged sterile packaging; dispose of it properly.
- Handle with care; store in a dry, dark environment below 25 °C.
- Do not re-sterilize any component.
- For use by healthcare professionals only.

PATIENT INFORMATION:

- Surgeons should inform patients about the limits of reconstruction and the need to protect the implant from full weight-bearing until adequate healing is achieved.
- Cement leakage during vertebroplasty may damage paravertebral structures, potentially causing spinal cord compression, intercostal neuralgia, cement leakage into the intervertebral space, perivertebral filling of veins and arteries (embolism risk), infection, and postoperative pain.
- To prevent cement leakage and detect undesired events in time, the procedure should be performed under real-time imaging guidance. Immediate operative intervention should be available to manage possible complications surgically.
- Preoperative comprehensive radiological evaluation is recommended to assess potential risks (e.g., vertebral body lesions, vertebral vascularization, or edema). Inadequate filling of the vertebral body with bone cement may lead to insufficient pain relief and compromise long-term stability of the treated vertebra.

PREPARATION AND APPLICATION:

After mixing the cement powder with the monomer liquid, a fast-setting, moldable dough is obtained. This dough is applied to bone cavities for stabilization and/or filling purposes.

Mixing the entire content of one bag of cement powder or container with the entire content of one ampoule of monomer liquid yields a single prepared dose. The required amount of cement dough depends on the specific surgical procedure and technique used.

Before starting the surgery, it is recommended to have at least one additional dose of bone cement prepared as a backup.

OGM VTB BONE CEMENT			
TEMPERATURE	20 °C	23 °C	26 °C
PREPARATION TIME	0' – 2'15"	0' – 2'	0' – 1'45"
OPERATING TIME	2'15" – 11'	2' – 10'30"	1'45" – 9'
HARDENING AND FINAL	11' – 17'	10'30" – 16'	9' – 14'

OPENING UNDER STERILE CONDITIONS:

The outer blister package must be opened under sterile conditions at the designated opening point, ensuring that the sterility of the inner powder package and glass ampoule is maintained. Before opening the powder package, gently shake the package or tap it lightly on a surface to allow the contents to settle at the bottom. This prevents powder loss. To facilitate opening of the glass ampoule, a predetermined breaking point is located at the neck between the ampoule body and head. To avoid glass fragments mixing with the cement, the ampoule should not be opened over the mixing container.

MIXING OF COMPONENTS:

The bone cement used during vertebroplasty can be applied to the vertebra using a percutaneous vertebroplasty or kyphoplasty system approved for controlled and stable injection. Refer to the manufacturer's instructions for the use of the system. Real-time X-ray monitoring (antero-lateral and lateral) is required during intravertebral application. If paravertebral cement leakage is detected, injection must be stopped immediately; injection may continue once the cement viscosity has sufficiently increased. If vertebral filling is insufficient, a separate contralateral entry can be applied.

After augmentation, the injection needle should be fitted with a mandrin, and after needle removal, any residual cement in the soft tissue should be removed.

Working and setting times depend on temperature, mixing speed, and environmental conditions. In particular, the temperatures of the cement powder, mixing system, table, and hands influence application times. As temperature increases, waiting, working, and setting times decrease.

STORAGE:

- Store the product in a dark, cool, and dry environment below 25 °C.
- Flammable; keep away from sources of ignition.
- Humidity should be between 0–50%.

SHELF LIFE / STERILITY:

The shelf life for sterile products is 36 months. The expiration date is printed on the folding box and protective blister (cap) package. Do not use the product if this date has passed. Opened or damaged protective blister (cap) packages or ampoule blisters cannot be re-sterilized and must be disposed of. Do not use the cement if the powder has turned yellow. The cement powder and packaging are sterilized by ethylene oxide, and the monomer liquid is sterilized by filtration.

PATIENT POPULATION:

Suitable for all patient populations except those with known or suspected hypersensitivity to the components of vertebroplasty cement, pregnant or breastfeeding patients, and patients with severe renal failure.

RESPONSIBILITY:

The manufacturer is not responsible for any direct, indirect, punitive, or related consequences arising from the reuse of vertebroplasty cement.

INTENDED USER:

Designed for use by healthcare professionals.

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